

Workshop on Observational Techniques

Basic introduction: spectrographs

Matti Dorsch¹


¹University of Potsdam

21. August 2023



at Ondřejov observatory




What do we want to achieve?



What do we want to achieve?



What do we want to achieve?



Compromise between


Large wavelength range

High efficiency

High spectral resolving power $R = \lambda/\Delta\lambda$

Accurate wavelength calibration

What do we want to achieve?




Compromise between

Large wavelength range

High efficiency

High spectral resolving power $R = \lambda/\Delta\lambda$

Accurate wavelength calibration



First-order spectroscopy

Long slit spectrograph






Figure: Schematic beam path in a long slit spectrograph.


A simple spectrograph




A simple spectrograph



A simple spectrograph




A simple spectrograph



Slits

Why do we need a slit?

Which slit width do we choose?




Slits

Why do we need a slit?

- take spectrum of *one* object

Which slit width do we choose?



Slits

Why do we need a slit?

- take spectrum of *one* object
- slit width $>$ PSF:
seeing-limited resolution
- slit width $<$ PSF:
slit-limited resolution
(also extended objects)
- Typically: 10 to 1000 μm

Which slit width do we choose?



Slits - light loss




Figure: Image scale for a lens with focal length f and infinite object distance.

$$\frac{\varphi}{2} \approx \tan \frac{\varphi}{2} = \frac{B}{2f} \quad \rightarrow \quad B = f \cdot \varphi$$

Typical seeing $\rightarrow \varphi \approx 2.5''$
 $f_{\text{Perek}} = 63.5 \text{ m}$ } Projected size on slit $B \approx 770 \mu\text{m}$

Slits - resolution

The slit used for OES is $b = 600 \mu\text{m} < 770 \mu\text{m}$! Why?





Figure: Schematic beam path in a long slit spectrograph.

Slits - resolution

The slit used for OES is $b = 600 \mu\text{m} < 770 \mu\text{m}$! Why?



Slits - resolution

The slit used for OES is $b = 600 \mu\text{m} < 770 \mu\text{m}$! Why?




Figure: Schematic beam path in a long slit spectrograph.

Reminder: $\Delta\alpha \stackrel{b \ll f_{\text{coll}}}{=} \frac{b}{f_{\text{coll}}}$, $n\lambda \stackrel{\text{interference}}{=} d \cdot (\sin \alpha + \sin \beta)$

Slits - resolution

The slit used for OES is $b = 600 \mu\text{m} < 770 \mu\text{m}$! Why?




Figure: Schematic beam path in a long slit spectrograph.

$$\text{Reminder: } \Delta\alpha \stackrel{b \ll f_{\text{coll}}}{\approx} \frac{b}{f_{\text{coll}}} \rightarrow \Delta\lambda \stackrel{\text{lin.}}{\approx} \frac{\partial\lambda}{\partial\alpha} \Delta\alpha = \frac{d}{n} \cos\alpha \Delta\alpha \stackrel{n=1}{\approx} \text{const.}$$

Slits - resolution

The slit used for OES is $b = 600 \mu\text{m} < 770 \mu\text{m}$! Why?





Figure: Schematic beam path in a long slit spectrograph.

$$\text{Reminder: } \Delta\alpha \stackrel{b \ll f_{\text{coll}}}{=} \frac{b}{f_{\text{coll}}} \rightarrow R_{\text{slit}} = \frac{\lambda}{\Delta\lambda_{\text{slit}}} = \frac{nf_{\text{coll}}}{db \cos \alpha} \lambda$$


Blaze grating - interference condition



$$n\lambda \stackrel{\text{interference}}{=} \Delta s = d \cdot (\sin \alpha + \sin \beta), \text{ order } n \in \mathbb{N}; \alpha + \beta \stackrel{\text{class.}}{=} 2\Theta_B$$

Δs : path difference between each groove

Blaze grating - diffraction orders



$$n\lambda_n^0 = d \cdot (\sin \alpha + \sin(2\Theta_B - \alpha)), \lambda_n^0 = \text{blaze wavelength (max. intensity)}$$

Blaze grating - blaze function





Figure: Dispersion and Blaze function $(\sin(x)/x)^2$ for a cross-dispersed échelle spectrograph.

Échelle spectrograph



- separate overlapping orders by cross-dispersion element
- optimized for high incidence angles and high orders:
 $\Theta_B = 69^\circ$ for OES

$$R_{\text{Échelle}} \approx \frac{f_{\text{coll}}}{b \cos \alpha} [\sin \alpha + \sin(2\Theta_B - \alpha)] \approx \text{constant} \approx 50000$$

Blaze grating - efficiency




Figure: Typical efficiency curves for blazed holographic gratings (edmundoptics).

CCD detector - efficiency




Figure: Quantum efficiency = % incident photons detected (SBIG ST-8XME).

Observation


Raw images




Raw images



Raw images



Raw images



Raw data summary

To produce a calibrated, 1-d spectrum, we need:

- *Science* frame

For each of these:

Raw data summary

To produce a calibrated, 1-d spectrum, we need:

- *Science* frame
- *Flat-field* frame

For each of these:

Raw data summary

To produce a calibrated, 1-d spectrum, we need:

- *Science* frame
- *Flat-field* frame
- *Calibration (arc)* frame

For each of these:

Raw data summary

To produce a calibrated, 1-d spectrum, we need:

- *Science* frame
- *Flat-field* frame
- *Calibration (arc)* frame

For each of these:

- *Bias* frame:
 - used to remove the CDD *readout signals*, including constant offset
 - taken with shortest exposure time and closed shutter
 - included in dark frame, required if $t_{\text{exp,dark}} \neq t_{\text{exp,science}}$

Raw data summary

To produce a calibrated, 1-d spectrum, we need:

- *Science* frame
- *Flat-field* frame
- *Calibration (arc)* frame

For each of these:

- *Bias* frame:
 - used to remove the CDD *readout signals*, including constant offset
 - taken with shortest exposure time and closed shutter
 - included in dark frame, required if $t_{\text{exp,dark}} \neq t_{\text{exp,science}}$
- *Dark* frame:
 - thermal excitation of electrons in the CCD leads to a constant background noise
 - also: hot/cold pixels/columns
 - taken with the same exposure time and temperature as science frame
 - has to be subtracted from science frame

Reduction steps

Bias frame




Figure: Median of 10 bias frames (closed shutter, shortest t_{exp} , log scale).

Bias frame




Figure: Median of 10 bias frames (closed shutter, shortest t_{exp} , log scale).

Dark frame / Cosmics




Figure: Dark frame detail, 3600s exposure (log scale).

Dark frame / Cosmics




Figure: Dark frame detail: median of five 3600s exposures (log scale).

Average frames




Figure: Science frame detail: 3600s exposure of BD+53 2790 (log scale).

Average frames




Figure: Science frame detail: dark frame subtracted.

Average frames




Figure: Science frame detail: dark frame subtracted, median of six exposures.

Sky background



Even the night sky is not completely black! Relevant for dark targets:

- air glow (emission lines due to chemical reactions in Earth's atmosphere, mainly at low altitudes $< 10^\circ$)
- scattered sunlight (astronomical twilight if Sun $< 18^\circ$ below horizon)
- moonlight
- *light pollution* (Potsdam, Berlin)
- in case of bad luck: planes (Tegel, Schönefeld)

Sky background




Figure: Potsdam sky background seen by DADOS (3h exposure average).

Sky background




Figure: Potsdam sky background seen by DADOS (3h exposure average).

Dispersion relation




Figure: NeAr calibration frame.

Find dispersion relation




Figure: Semi-automatic emission line identification (NeXe lamp).

Find dispersion relation




Figure: Dispersion relation as deviation from a linear relation between pixel and wavelength (NeXe lamp).

Find dispersion relation




Figure: Deviations from the fitted dispersion relation (NeXe lamp).

Step-by-step summary




Figure: DADOS spectrum of Alcyone: dark, averaged.

Step-by-step summary




Figure: DADOS spectrum of Alcyone: dark, averaged, flat.

Step-by-step summary




Figure: DADOS spectrum of Alcyone: dark, averaged, flat, calibrated.

Optional: (relative) flux calibration

Spectral shape is still affected by the flatfield shape, etc ...

+

We still don't have physical units

Optional: (relative) flux calibration

To fix this, create a standard star:

- take spectrum of a laboratory source with known flux distribution $L(\lambda)$:

$$R_{\text{lab}}(\lambda)$$

- then, the calibration factor is:

$$C(\lambda) = L(\lambda)/R_{\text{lab}}(\lambda)$$

- take spectrum of a standard star (correct for atm. extinction!):

$$R_{\text{std}}(\lambda)$$

- then, the flux of the std. star in physical units ($\text{erg cm}^{-2} \text{s}^{-1} \text{\AA}^{-1}$) is:

$$f_{\text{std}}(\lambda) = C(\lambda)R_{\text{std}}(\lambda) = \frac{R_{\text{std}}(\lambda)}{R_{\text{lab}}(\lambda)}L(\lambda)$$

Optional: (relative) flux calibration

If flux distribution of a standard (comparison) star $f_{\text{std}}(\lambda)$ is known:

- take high S/N spectrum of standard star:

$$R_{\text{std}}(\lambda)$$

- then, the "calibration" factor is:

$$c(\lambda) = f_{\text{std}}(\lambda)/R_{\text{std}}(\lambda)$$

- take spectrum of target star $R_*(\lambda)$, then:


$$f_*(\lambda) = c(\lambda)R_*(\lambda) = \frac{R_*(\lambda)}{R_{\text{std}}(\lambda)}f_{\text{std}}(\lambda)$$

- the best calibrated star is Vega (no atm. extinction):

$$f_{\text{Vega}}(5556 \text{ \AA}) = 3.44 \pm 0.05 \text{ erg cm}^{-2} \text{ s}^{-1} \text{ \AA}^{-1}$$


Optional: (relative) flux calibration

Here, relative flux calibration using excellent synthetic spectra:




Optional: (relative) flux calibration

Here, relative flux calibration using excellent synthetic spectra:




Optional: (relative) flux calibration

Here, relative flux calibration using excellent synthetic spectra:




Optional: (relative) flux calibration

Here, relative flux calibration using excellent synthetic spectra:




Optional: (relative) flux calibration

Here, relative flux calibration using excellent synthetic spectra:




Optional: (relative) flux calibration

Here, relative flux calibration using excellent synthetic spectra:




Optional: (relative) flux calibration

Here, relative flux calibration using excellent synthetic spectra:




Optional: (relative) flux calibration

Here, relative flux calibration using excellent synthetic spectra:



Optional: (relative) flux calibration

Here, relative flux calibration using excellent synthetic spectra:



Why not absolute flux calibration?

Absolute flux calibration would be even better!

Fit synthetic spectra to get:


- angular diameter (solid angle) $\Theta = 2R/D$
- temperature T_{eff}
- color excess (interstellar reddening) E_{B-V}
- using parallax $\varpi \rightarrow$ stellar radius $R = \Theta/(2\varpi)$

Problems:

- airmass changes with distance to horizon
→ wavelength dependent extinction
- star may leave the slit during the observation
- clouds

Optional: determine resolution with calibration lamp

We need to know the spectral resolution to fit spectra!



Arc lines are intrinsically sharp \rightarrow fit Gaussian: $R = \lambda_{\text{line}}/\text{FWHM}_{\text{line}}$.

Example reduction using python scripts

In three steps (adjust file names, run scripts with -h for help):

```
# stack hgar images and flat images for star
~/scripts/evolved/0_average_images.py  hgar/ -o
    star_hgar_stacked.fit
```

```
# identify calibration lines
~/scripts/evolved/1_findcaliblines.py  -arc
    star_hgar_stacked.fit -rsc 500 570
```

```
# apply calibration and extract spectra
~/scripts/evolved/2_extractspectrum.py -sc
    star_1200s_stacked.fit -df ~/data/20190903/dark
    /1200s/ -ff flats/ -fd ~/data/20190903/dark/1s/
    -rsc 500 570
```

Questions?